

Supplemental Table S1. Antibody Information

Antibody (code)	Company	Staining kit	Autostainer	Dilution
Gal-3 (NCL-GAL3)	Novo	OptiView DAB Detection Kit	BenchMark XT	1:200
HBME-1 (M3505)	Dako	OptiView DAB Detection Kit	BenchMark XT	1:200
TTF-1 (NCL-TTF-1)	Novo	OptiView DAB Detection Kit	BenchMark XT	1:200
Mitochondria antibody (AB92824)	ABCAM	OptiView DAB Detection Kit	BenchMark XT	1:2,000
Thyroglobulin (M0781)	Dako	OptiView DAB Detection Kit	BenchMark XT	1:4,000
PAX-8 (363M-15)	MRQ50	OptiView DAB Detection Kit	BenchMark XT	1:100

Gal-3, galectin-3; HBME-1, Hector Battifora mesothelial cell-1; TTF-1, thyroid transcription factor-1; PAX-8, paired box gene 8.